



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Timer,60 min

**Product Code :**  
LBNY-0018-2200012



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Countdown timer.  
Plastic outer casing.  
Robust design.  
Loud sound to indicating set time has elapsed.  
Easily readable dial with large numbers.  
Smooth surface which allows for easy cleaning.  
Mechanical, spring driven.  
Range: 0 - 60 minutes.  
Size allows for easy transporting.

**Technical Specification :**

WEIGHT AND VOLUME (Packaged)  
Volume: 0.300 dm<sup>3</sup>.  
Weight: 77 gr.



**Laboratory Instrument India**

