



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Papersheet, crepe, for ster. pack

**Product Code :**

LBNY-0018-2200014



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

It is designed to allow sterilization of the enclosed medical device and also to maintain sterility of the device until the packaging is opened for use of the device, or until a predetermined shelf date is expired. This is a single-use device. Paper sheets intended to be used to contain medical devices that are to be sterilized.

**Technical Specification :**

Good permeability to steam.  
Good microbiological barrier.  
Low permeability to air.  
Water repellence time: at least 20 seconds.  
Non-sterile and single use.  
Cellulose based.  
Grammage: 58 - 62 g/m<sup>2</sup>  
Porosity: 50µm maximum.  
Size of the individual sheets: 0.9 x 0.9 m



**Laboratory Instrument India**

