



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Retractor, abdo, Balfour, 3 blades

**Product Code :**

LBNY-0019-220001



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Surgical instrument to separate fatty tissue, muscles or viscera after the abdominal incision to expose the operative field for examination or intervention.

Static function instrument: abdominal self-retaining retractor.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.50% carbon, 12-15% chromium, 0.45-0.90% molybdenum.

The compass (frame) is a support with hinged arms.

2 fenestrated end blades: approx. 35 x 70 mm.

1 center blade, adjustable in height: approx. 45 x 80 mm.

See Martin Catalogue no. 24 Ref: 15 852 00.

Locking mechanism and compass (self-holding).

Maximum wound opening: approx. 180 mm.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 40% RH to 70% RH.

Storage and transport conditions: 5°C to 40°C.

---

# Laboratory Instrument India

