



Email : sales@laboratoryinstrumentindia.com

Product Name :

Clamp,towel,Backhaus,130mm

Product Code :

LBNY-0019-220002



www.laboratoryinstrumentindia.com

Description :

Pressure force instrument: towel clamp, springy.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.
One fine single tooth at each jaw.
Reusable.
Hard ratchet, lockable.
Length: 13 cm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India