



Email : sales@laboratoryinstrumentindia.com

Product Name :

Curette,uterine,Sims,8mm,blunt

Product Code :

LBNY-0019-220003



www.laboratoryinstrumentindia.com

Description :

Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.16 - 0.20 % carbon; 13 % chromium.
Blade shape: Oval, fenestrated.
Blade edge type: Blunt.
Uterine curette.
Shank type: Rigid.
Length: shaft 26 cm.
Width: spoon 8 mm.
Reusable

Technical Specification :

Environmental conditions:
Humidity: 40 % - 70 % RH
Storage and transport conditions: 5°C - 40°C



Laboratory Instrument India

