



Email : sales@laboratoryinstrumentindia.com

Product Name :

Curette,uterine,Sims,7mm,sharp

Product Code :

LBNY-0019-220004



www.laboratoryinstrumentindia.com

Description :

For removal of retained products of conception from inside the uterus.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15 - 0.40 % carbon; 13 - 14 % chromium.
Uterine curette.
Blade edge type: Sharp.
Shank type: Rigid.
Blade shape: Oval, fenestrated.
Length: shaft approx. 26 cm.
Width: spoon 7 mm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 40 % RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

