



Email : sales@laboratoryinstrumentindia.com

Product Name :

Curette,uterine,Sims,9mm,sharp

Product Code :

LBNY-0019-220005



www.laboratoryinstrumentindia.com

Description :

Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.16-0.40 % carbon; 13-14 % chromium.
Blade shape: Oval, fenestrated.
Blade edge type: Sharp.
Uterine curette.
Shank type: Rigid.
Reusable.
Length: Shaft 26 cm.
Width: Spoon 9 mm.

Technical Specification :

Environmental conditions:
Humidity at 40% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

