



Email : sales@laboratoryinstrumentindia.com

Product Name :

Curette,uterine,Sims,12mm,sharp

Product Code :

LBNY-0019-220006



www.laboratoryinstrumentindia.com

Description :

For removal of retained products of conception from inside the uterus.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.40 % carbon, 13-14% chromium.
Uterine curette.
Blade edge type: Sharp.
Length: Shaft 26 cm.
Width: Spoon 12 mm.
Shank type: Rigid.
Blade shape: Oval, fenestrated.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 40% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

