



Email : sales@laboratoryinstrumentindia.com

Product Name :

Scoop,uterine,Simon,12mm,sharp

Product Code :

LBNY-0019-220007



www.laboratoryinstrumentindia.com

Description :

For removal of retained products of conception from inside the uterus during gynaecological procedures.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.16-0.40 % carbon; 14 % chromium.

Shank type: Rigid.

Blade shape: Oval form of a spoon.

Handle: Austenitic steel

Blade edge type: Sharp

Length: Shaft 25 cm.

Width: Spoon 12 mm.

Reusable.

Technical Specification :

Environmental conditions:

Humidity: 40% RH to 70% RH.

Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

