



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,dressing,standard,155mm,str

Product Code :

LBNY-0019-220009



www.laboratoryinstrumentindia.com

Description :

Used in surgery and nursing activities.

Forceps without teeth are used for dissecting delicate tissues.

Dressing forceps for gripping, dissecting tissue and coagulation of vessels.

Technical Specification :

Pressure force instrument: dissecting forceps, springy

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.20-0.40% carbon; 13-14% chromium.

Straight without teeth.

Serrated jaws

Flexible arms.

Length: 15.5 cm.

Reusable.



Laboratory Instrument India

