



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps,dressing,standard,155mm,str

**Product Code :**

LBNY-0019-220009



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Used in surgery and nursing activities.

Forceps without teeth are used for dissecting delicate tissues.

Dressing forceps for gripping, dissecting tissue and coagulation of vessels.

**Technical Specification :**

Pressure force instrument: dissecting forceps, springy

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.20-0.40% carbon; 13-14% chromium.

Straight without teeth.

Serrated jaws

Flexible arms.

Length: 15.5 cm.

Reusable.



**Laboratory Instrument India**

