



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,dressing,standard,250mm,str

Product Code :

LBNY-0019-2200010



www.laboratoryinstrumentindia.com

Description :

Dressing forceps for gripping, dissecting tissue and coagulation of vessels.
Forceps without teeth are used for dissecting delicate tissues.
Used in surgery and nursing activities.

Technical Specification :

Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.20% carbon, 13% chromium.
Straight without teeth.
Serrated jaws
Flexible arms.
Length: 25 cm.
Reusable.



Laboratory Instrument India