



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,dressing,Cheron,250mm

Product Code :

LBNY-0019-2200011



www.laboratoryinstrumentindia.com

Description :

Forceps for dressing / swabbing during gynaecological procedures.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.
Reusable.
Variable setting of ratchet, lockable.
Adjustment of jaws.
Scissors-like hinged design with non-cutting blades.
Flexible arms.
Length: 25 cm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

