




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Flat Files Of Different Lengths

**Product Code :**  
LBNY-0004-820028



**Description :**

Flat Files Of Different Lengths

**Technical Specification :**

Flat Files Of Different Lengths  
20 Pcs Metal File Set  
1ea \* 4 Way wood rasp file  
1ea \* Flat file(8"/200mm)  
1ea \* Half-round file(8"/200mm)  
1ea \* Round file(8"/200mm)  
1ea \* Triangle File(8"/200mm)  
14ea \* Needle file (4,2"/108mm)  
1ea \* Cleaning Brush  
1ea x Portable Carryhtg bag  
All the 20pcs file set are made of high carbon forged steel material T12.  
Material : Stainless steel

---

# Laboratory Instrument India

