



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,artery,Kelly,140mm,cvd

Product Code :

LBNY-0019-2200012



www.laboratoryinstrumentindia.com

Description :

For haemostasis, clamping of large blood vessels or manipulating heavy tissue.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.

Reusable.

Variable settings of the ratchet, lockable.

Adjustment of the jaws.

Pressure force instrument: haemostatic forceps, springy.

Curved.

Flexible arms.

Length: 14 cm.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

