



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps,artery,Kelly,140mm,cvd

**Product Code :**

LBNY-0019-2200012



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

For haemostasis, clamping of large blood vessels or manipulating heavy tissue.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.

Reusable.

Variable settings of the ratchet, lockable.

Adjustment of the jaws.

Pressure force instrument: haemostatic forceps, springy.

Curved.

Flexible arms.

Length: 14 cm.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

