



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,artery,Mixer,230mm

Product Code :

LBNY-0019-2200013



www.laboratoryinstrumentindia.com

Description :

For dissecting, clamping, grasping tissues and for passing a thread around a vein or artery.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.20% carbon; 13% chromium.

Reusable.

Flexible ratchet, lockable.

Jaws which grip the thread well.

Pressure force instrument: threading forceps, springy.

Length: 23 cm.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

