



Email : sales@laboratoryinstrumentindia.com

Product Name :
Forceps,artery,Mixer,230mm

Product Code :
LBNY-0019-2200013



www.laboratoryinstrumentindia.com

Description :

For dissecting, clamping, grasping tissues and for passing a thread around a vein or artery.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.20% carbon; 13% chromium.
Reusable.
Flexible ratchet, lockable.
Jaws which grip the thread well.
Pressure force instrument: threading forceps, springy.
Length: 23 cm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

