



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps,artery,Kocher,140mm,str

**Product Code :**

LBNY-0019-2200014



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Multi-purpose use in surgery: For haemostasis, gripping, dissection

Variable setting of the ratchet, lockable.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.20% carbon, 13% chromium.

Adjustment of the jaws.

Pressure force instrument: haemostatic forceps, springy.

Straight.

Flexible arms.

Reusable.

Length: 14 cm.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 5°C to 40°C.



**Laboratory Instrument India**

