



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,artery,Kocher,140mm,str

Product Code :

LBNY-0019-2200014



www.laboratoryinstrumentindia.com

Description :

Multi-purpose use in surgery: For haemostasis, gripping, dissection
Variable setting of the ratchet, lockable.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.20% carbon, 13% chromium.
Adjustment of the jaws.
Pressure force instrument: haemostatic forceps, springy.
Straight.
Flexible arms.
Reusable.
Length: 14 cm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

