



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,artery,Pean/Roch,200mm,cvd

Product Code :

LBNY-0019-2200015



www.laboratoryinstrumentindia.com

Description :

For clamping larger tissue and vessels for haemostasis in multiple surgical procedures.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.

Adjustment of jaws.

Length: approx. 20 cm.

Reusable.

Pressure force instrument, haemostatic forceps, springy.

Curved.

Variable setting of ratchet, lockable.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

