



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,artery,Pean/Roch,240mm,cvd

Product Code :

LBNY-0019-2200016



www.laboratoryinstrumentindia.com

Description :

Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.
For clamping larger tissue and vessels for haemostasis in multiple surgical procedures.
Pressure force instrument: haemostatic forceps, springy.
Variable setting of ratchet, lockable.
Adjustment of jaws.
Curved.
Length: approx. 24 cm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

