



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,intest,clamp,Doyen,230mm,cvd

Product Code :

LBNY-0019-2200017



www.laboratoryinstrumentindia.com

Description :

For transverse occlusion of a section of intestine while suturing, resection or anastomosis.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.

Reusable.

Occlusion: The form and mechanical action must permit precise closing.

Adhesion: Provided by the surface of the jaws.

Smooth, springy jaws with atraumatic longitudinal serrations.

Length: 23 cm.

Pressure force instrument: intestinal clamp, springy

Non-crushing clamp, curved.

Long jaws and shorter shanks.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

Laboratory Instrument India

