



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps,artery,Halst-Mosq,125mm,cvd

**Product Code :**

LBNY-0019-2200018



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

For clamping vessels for haemostasis in general surgical procedures.  
Material: Martensitic steel (quenched, magnetic steel).  
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.  
Pressure force instrument: haemostatic forceps, springy.  
Slim, curved forceps without teeth.  
Flexible arms.  
Atraumatic jaws.  
Variable setting of the ratchet, lockable.  
Length: 12.5 cm.  
Reusable.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

