



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,artery,Halst-Mosq,125mm,cvd

Product Code :

LBNY-0019-2200018



www.laboratoryinstrumentindia.com

Description :

For clamping vessels for haemostasis in general surgical procedures.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.
Pressure force instrument: haemostatic forceps, springy.
Slim, curved forceps without teeth.
Flexible arms.
Atraumatic jaws.
Variable setting of the ratchet, lockable.
Length: 12.5 cm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

