



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,tissue,Babcock,200mm

Product Code :

LBNY-0019-2200019



www.laboratoryinstrumentindia.com

Description :

To grip soft and delicate tissues (particularly intestines)
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.
Reusable.
Soft ratchet, lockable.
Pronounced but atraumatic ridges of the grippers.
Ridged grippers with aperture (triangular).
Pressure force instrument: grasping forceps, springy
Flexible arms.
Length: 20 cm.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

