



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,tissue,standard,145mm,str

Product Code :

LBNY-0019-2200020



www.laboratoryinstrumentindia.com

Description :

Used in surgery for gripping, dissecting thick tissues and coagulation of vessels.

Adjustment and good gripping of the jaws.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.

Reusable.

Pressure force instrument: dissecting forceps, springy.

Straight.

With 1 x 2 teeth.

Flexible arms.

Length: 14.5 cm.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

