



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,tissue,Collin,160mm

Product Code :

LBNY-0019-2200021



www.laboratoryinstrumentindia.com

Description :

For gripping soft and delicate tissues (lungs and intestines).
Material: Martensitic steel (quenched, magnetic steel)
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.
Pressure force instrument: grasping forceps, springy.
Atraumatic ridged grippers, each with aperture (oval).
Soft ratchet, lockable.
Length: 16 cm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

