



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Pipe Bender

Product Code :
LBNY-0004-820029



Description :

Hydraulic Pipe Bender

Technical Specification :

Hydraulic Pipe Bender

Suitable for bending of low pressure fluid transport welding pipe, arc welded steel pipe, ordinary steel pipe under cold condition. It is not suitable for thin wall pipe and non-ferrous metal pipe

Before use, it is necessary to loosen the gold bronze oil valve of the cylinder so as not to break the seal ring.

Before use, it is necessary to screw a black screw on the cylinder so that air in the cylinder can be discharged.

Ductile cast iron dies, strength, toughness and plasticity are good, durable and not easily broken

Working Pressure PSI 16ton

Recommended Dimension : L x W x H n. 24 3/4 x 8 3/4 x 22 1/2

Top rollers adjust for 90 bend,

Bending Applications Schedule 40, 80 pipe

Six Cast Iron Shoes



Laboratory Instrument India

