



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps,tissue,standard,250mm,str

**Product Code :**

LBNY-0019-2200022



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Used in surgery for gripping, dissecting thick tissues and coagulation of vessels.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.

Adjustment and good gripping of the jaws.

Pressure force instrument: dissecting forceps, springy.

Straight.

With 1 x 2 teeth.

Flexible arms.

Length: 25 cm.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

