



Email : sales@laboratoryinstrumentindia.com

Product Name :
Forceps,tissue,Allis,150mm

Product Code :
LBNY-0019-2200023



www.laboratoryinstrumentindia.com

Description :

For grasping organs, soft or heavy tissue.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon, 12-13% chromium.
Precise adjustment of the teeth.
Hard ratchet, lockable.
Pressure force instrument: tissue grasping forceps, springy.
4 x 5 teeth.
Length: 15 cm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

