



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Forceps,tissue,Duval,230mm

**Product Code :**  
LBNY-0019-2200024



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Martensitic steel (quenched, magnetic steel).  
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.  
To grip soft and delicate tissues (lungs, intestines, stomach).  
Pronounced but atraumatic ridges of the grippers.  
Ridged grippers with aperture (triangular).  
Pressure force instrument: gripping forceps, springy  
Soft ratchet, lockable.  
Flexible arms.  
Length: 23 cm.  
Jaw: approx. 27 mm.  
Reusable.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.

---

# Laboratory Instrument India

