



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps,tissue,Duval,230mm

Product Code :

LBNY-0019-2200024



www.laboratoryinstrumentindia.com

Description :

Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20% carbon; 12-13% chromium.
To grip soft and delicate tissues (lungs, intestines, stomach).
Pronounced but atraumatic ridges of the grippers.
Ridged grippers with aperture (triangular).
Pressure force instrument: gripping forceps, springy
Soft ratchet, lockable.
Flexible arms.
Length: 23 cm.
Jaw: approx. 27 mm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.

Laboratory Instrument India

