



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps, uterine, Duplay, 280mm, cvd

**Product Code :**

LBNY-0019-220025



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

For gripping and immobilisation of thick tissue (uterine traction forceps).

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0-15-0.20% carbon, 12-13% chromium.

Hard ratchet, lockable.

Pressure force instrument: gripping forceps, springy.

Uterine traction forceps, vulsellum / tenaculum, double curved.

Length: 28 cm.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

