



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps, uterine, Duplay, 280mm, cvd

Product Code :

LBNY-0019-2200025



www.laboratoryinstrumentindia.com

Description :

For gripping and immobilisation of thick tissue (uterine traction forceps).

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0-15-0.20% carbon, 12-13% chromium.

Hard ratchet, lockable.

Pressure force instrument: gripping forceps, springy.

Uterine traction forceps, vulsellum / tenaculum, double curved.

Length: 28 cm.

Reusable.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

