



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps, uterine, Phaneuf, 215mm, cvd

Product Code :

LBNY-0019-2200026



www.laboratoryinstrumentindia.com

Description :

For haemostasis of the arteries, especially uterine arteries.
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.20 % carbon, 12-13 % chromium.
Pressure force instrument: Hemostatic forceps, slightly springy.
Adjustment of the jaws.
Curved, angled on flat forceps with 1 x 2 teeth.
Flexible arms.
Variable setting of ratchet (lockable).
Length: 21.5 cm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 30% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

