



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps, uterine, Phaneuf, 215mm, cvd

**Product Code :**

LBNY-0019-2200026



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

For haemostasis of the arteries, especially uterine arteries.  
Material: Martensitic steel (quenched, magnetic steel).  
Martensitic steel composition: 0.15-0.20 % carbon, 12-13 % chromium.  
Pressure force instrument: Hemostatic forceps, slightly springy.  
Adjustment of the jaws.  
Curved, angled on flat forceps with 1 x 2 teeth.  
Flexible arms.  
Variable setting of ratchet (lockable).  
Length: 21.5 cm.  
Reusable.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

