



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Forceps, uterine, Museux, 240mm, cvd

**Product Code :**

LBNY-0019-2200027



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Used to grip and to immobilise thick tissue (uterine traction forceps) during obstetrical/gynaecological procedures.

Precise adjustment of the teeth.

Hard ratchet, lockable.

Pressure force instrument: gripping forceps, springy

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0-15-0.20 % carbon; 12-13 % chromium.

Uterine traction forceps, vulsellum / tenaculum, curved sideways.

With 2 x 2 teeth.

Length: 24 cm.

Jaws: 8-10 mm.

Reusable.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

---

# Laboratory Instrument India

