



Email : sales@laboratoryinstrumentindia.com

Product Name :

Forceps, uterine, Museux, 240mm, cvd

Product Code :

LBNY-0019-2200027



www.laboratoryinstrumentindia.com

Description :

Used to grip and to immobilise thick tissue (uterine traction forceps) during obstetrical/gynaecological procedures.

Precise adjustment of the teeth.

Hard ratchet, lockable.

Pressure force instrument: gripping forceps, springy

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0-15-0.20 % carbon; 12-13 % chromium.

Uterine traction forceps, vulsellum / tenaculum, curved sideways.

With 2 x 2 teeth.

Length: 24 cm.

Jaws: 8-10 mm.

Reusable.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.

Laboratory Instrument India

