



Email : sales@laboratoryinstrumentindia.com

Product Name :

Needle holder, Mayo-Hegar, 180mm, str

Product Code :

LBNY-0019-2200028



www.laboratoryinstrumentindia.com

Description :

Surgical instrument to hold a suture needle while stitching through tissue.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20 % carbon, 12-13 % chromium.

A well-defined longitudinal groove to prevent deformation of the needle.

Jaws with pronounced cross ridges.

Pressure force instrument: needle holder, springy.

Straight with ring handles.

Progressive setting of the ratchet to enable the needle to be gripped with varying tightness.

Reusable.

Length: 18 cm.

Technical Specification :

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

