



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Needle holder, Mayo-Hegar, 180mm, str

**Product Code :**

LBNY-0019-2200028



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Surgical instrument to hold a suture needle while stitching through tissue.

Material: Martensitic steel (quenched, magnetic steel).

Martensitic steel composition: 0.15-0.20 % carbon, 12-13 % chromium.

A well-defined longitudinal groove to prevent deformation of the needle.

Jaws with pronounced cross ridges.

Pressure force instrument: needle holder, springy.

Straight with ring handles.

Progressive setting of the ratchet to enable the needle to be gripped with varying tightness.

Reusable.

Length: 18 cm.

**Technical Specification :**

Environmental conditions:

Humidity: 30% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

