



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hand Pipe Bending Tool

Product Code :
LBNY-0004-820030



Description :

Hand Pipe Bending Tool

Technical Specification :

Hand Pipe Bending Tool

High Quality: Made of high quality carbon steel, heavy quality manual pipe bender for pipe bending, it has a long service life

Easy Storage: Packed in heavy-duty wooden cases for better protection and easy storage, that can avoid tool damage

Easy to Use: It is easy to operate and it is suitable for bending steel tubes, brass, copper, wrought iron and aluminum tubes.

Included Components:

Bender body,

Crossbar

Reverse bending

Attachment and mandrels for 1/4", 5/16", 3/8", 7/2", 5/8", 3/4" and 7"8 tubing

Pump lever until angle up to 90 degrees

Laboratory Instrument India

