



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Probe, double-ended, 145mm

**Product Code :**  
LBNY-0019-2200034



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Austenitic steel (non-quenched, non-magnetic steel).  
Austenitic steel composition: 16-18% chromium, 8-14% nickel, 2-3% molybdenum.  
To explore a surgical wound or to explore and/or dilate strictures inside a fistula.  
Probe / dilator, straight.  
Blunt bulbous shaped ends.  
Diameter of the ends: approx. 2 mm.  
Length: approx. 145 mm.  
Reusable.

**Technical Specification :**

Environmental conditions:  
Humidity: 30% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

