



Email : sales@laboratoryinstrumentindia.com

Product Name :
Tube suction,270mm

Product Code :
LBNY-0019-2200035



www.laboratoryinstrumentindia.com

Description :

Used for aspirating fluids / secretions from the operative area.
Material: Austenitic steel (non-quenched, non-magnetic steel).
Austenitic steel composition: 16-18% chromium, 8-14% nickel, 2-3% molybdenum.
Curved.
Suction tube with removable tip and tubing connector.
Length: approx. 27 cm.
Diameter: approx. 6 mm or 8 mm.
Reusable.

Technical Specification :

Environmental conditions:
Humidity: 40% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

