



Email : sales@laboratoryinstrumentindia.com

Product Name :
Retractor, abdo, Collin, 3 blades

Product Code :
LBNY-0019-2200036



www.laboratoryinstrumentindia.com

Description :

To retract skin, fatty tissue, muscles or viscera after the abdominal incision to expose the operative field.
3 blades, the two lateral blades being mobile and the center blade removable.
Static function instrument: abdominal self-retaining retractor.
With locking mechanism (self-holding).
Material: Martensitic steel (quenched, magnetic steel).
Martensitic steel composition: 0.15-0.50% carbon, 12-15% chromium, 0.45-0.90% molybdenum.
Reusable.
Solid lateral blades: approx. 38 x 60 mm.
Solid center blade: approx. 38 x 55 mm.
Maximum wound opening: approx. 100 mm.

Technical Specification :

Environmental conditions:
Humidity: 40% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

