



Email : sales@laboratoryinstrumentindia.com

Product Name :
Retractor,vaginal,45x85mm

Product Code :
LBNY-0019-2200037



www.laboratoryinstrumentindia.com

Description :

Static function instrument: vaginal retractor with handle.
Material: Austenitic steel (non-quenched, non-magnetic steel).
Austenitic steel composition: 18 to 20% chromium, 8 to 10% nickel.
Reusable.
Blunt lateral edges.
Blade length: approx. 85 mm.
Blade width: approx. 45 mm.
Length: approx. 24 cm.

Technical Specification :

Environmental conditions:
Humidity: 40% RH to 70% RH.
Storage and transport conditions: 5°C to 40°C.



Laboratory Instrument India

