



Email : sales@laboratoryinstrumentindia.com

Product Name :
Potable Tool Set

Product Code :
LBNY-0004-820031



Description :

Potable Tool Set (Total: 56 pcs)

Technical Specification :

Potable Tool Set (Total: 56 pcs)
Socket(6pt) : 45-08, 10, 12,13,14,17,19,21,22,24
Plug Socket: 4P-16,20.6
Extension Bar : 320(15m),321(150mm)
Long spanner: LS-5.507, 0810, 1012, 1113,1214,1417
Mechanical Offset Wrench Long Type: M44-0810, 1012, 7417, 7921, 2224
Short Offset Wrench:M46-1012
Flare Nut Wrench: Nn6-1012
Adjustable Wrench : MW-300
Needle Nose Pliers; RP-150
Combination pliers: CP'200
Cutting Pliers: KN-150
Power screwdriver: PSD-M2, PSD-P2
Heavy Duty Screwdriver: YMD-150
Ball Pin flammer: BH-10 (1 pound)
Brake Shu Adjusting Tool : BST:-1 10
Chisel: 19X165mm
Scrapct:25X160mm
Punch: 125mm

Ball Point Hex Key Wrench Set; BPL900 (1.5, 2, 2.5, 3, 4, 5, 6, 8, 10)
Exclusive tool case (2215)



Laboratory Instrument India