



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sound,uterine,Martin,320mm

Product Code :
LBNY-0019-2200049



www.laboratoryinstrumentindia.com

Description :

Material: Silver-coated copper with brass handle or Austenitic steel (non-quenched, non-magnetic steel).
Austenitic steel composition: 16 - 18% chromium; 10 - 14% nickel; 2 - 3% molybdenum.
To measure and/or explore the depth or length of the uterus, cervix, and vagina.
Bulbous distal end, mounted on a handle.
Static function instrument: uterine sound.
Malleable shaft graduated in cm.
Reusable.
Length: approx. 32 cm.

Technical Specification :

Environmental conditions:
Humidity: 40% RH to 70% RH.
Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

