



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sound,uterine,Martin,320mm

**Product Code :**  
LBNY-0019-2200049



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Material: Silver-coated copper with brass handle or Austenitic steel (non-quenched, non-magnetic steel).  
Austenitic steel composition: 16 - 18% chromium; 10 - 14% nickel; 2 - 3% molybdenum.  
To measure and/or explore the depth or length of the uterus, cervix, and vagina.  
Bulbous distal end, mounted on a handle.  
Static function instrument: uterine sound.  
Malleable shaft graduated in cm.  
Reusable.  
Length: approx. 32 cm.

**Technical Specification :**

Environmental conditions:  
Humidity: 40% RH to 70% RH.  
Storage and transport conditions: 1°C to 30°C.



**Laboratory Instrument India**

