



Email : sales@laboratoryinstrumentindia.com

Product Name :

Speculum,vaginal,Graves,75x20mm

Product Code :

LBNY-0019-220050



www.laboratoryinstrumentindia.com

Description :

Material: Austenitic steel (non-quenched, non-magnetic steel).

To dilate the vagina to examine the vaginal canal and cervix.

Speculum, vaginal, self-retaining.

Double blade with a screw mechanism to keep the blades open.

Austenitic steel composition: 16 - 18% chromium; 10 - 14% nickel; 2 - 3% molybdenum.

Blade length: approx. 75 mm.

Blade width: approx. 20 mm.

Reusable.

Technical Specification :

Environmental conditions:

Humidity: 40% RH to 70% RH.

Storage and transport conditions: 1°C to 30°C.



Laboratory Instrument India

